



Technical Data Sheet

OSSFILATM PEEK

OSF-FLC

3D Printer Filament

OssfilaTM **PEEK** is a medical grade 3D printing material suitable for creating patient specific orthopaedic implants or surgical tools.

The production of Ossfila[™] PEEK follows ISO 13485 and ISO 14644. The raw materials meet the requirements of ISO 10993. The final products meet the requirements of USP Class VI and ISO 10993-6:2016.

PRODUCT FEATURES

- USP Class VI certified
- ISO 10993-6:2016 certified
- Excellent biocompatibility
- Excellent mechanical strength
- Non-toxic
- Stable to autoclave, gamma radiation,
 E-beam, and ethylene oxide

APPLICATIONS

OssfilaTM PEEK is recommended for creating patient specific orthopaedic implants, spinal implants, dental implants, and surgical tools.

IMPORTANT NOTICE

- The Filament is intended to be used as a raw material to additively manufacture medical devices, surgical tools and/or equipment parts via an FDM 3D printer.
- Store the Filament in a dry cool place and away from heat. Absorption of moisture could affect filament quality.
- All recommendations in this document are based on the operations conducted by the manufacturer's
 facility. While every effort has been made to ensure safety, accuracy and reliability, it is imperative
 that the user verifies and/or validates the process parameters with their own facility's equipment and
 conditions if deemed necessary.





PRODUCT SPECIFICATIONS				
Key Materials	Polyether ether ketone			
Filament Diameter	1.75 mm ± 0.05 mm			
Net Weight	500 g			
Colour	Natural			
Storage Condition	Sealed, cool and dry place			
Drying Condition	120°C for 4 hours			
Features	Biocompatible, medical grade			

RECOMMENDED PRINTER SETTINGS				
Extrusion/ Nozzle	390 – 410°C			
Temperature				
Build Plate	130 – 150°C			
Temperature	130 – 130 C			
Chamber	90°C or higher as appropriate			
Temperature				
Adhesion Type	Glue or Kapton® film			
Nozzle Diameter	> 0.4 mm			
	_			

MECHANICAL PROPERTIES					
	XZ Orientation	XY Orientation	Method		
Tensile Properties					
Ultimate Tensile Strength	78 MPa	86 MPa	ASTM D638 (Type IV)		
Young's Modulus	3030 MPa	3250 MPa	ASTM D638 (Type IV)		
Elongation	6%	5%	ASTM D638 (Type IV)		
Flexural Properties					
Flexural Strength	146 MPa	151 MPa	ASTM D790		
Flexural Modulus	3110 MPa	3570 MPa	ASTM D790		

THERMAL PROPERTIES					
	Test Value	Method			
Melting Point	343°C	ISO 11357			
Glass Transition Temperature	143°C	ISO 11357			





Note:

- 1. The test articles of the mechanical test are 3D printed objects printed to the dimension specified in ASTM D638 and ASTM D790.
- 2. Test articles were printed with 100% infill and line as the infill pattern.
- 3. It is recommended that the 3D printer equips with a closed chamber and a heating build plate which could reach at least 130°C.
- 4. It is recommended to use a raft instead of a brim for build plate adhesion.
- 5. It is recommended to dry the Filament at 120°C for 4 hours before printing for best performance. If the mentioned conditions cannot be achieved, dry the filament at 50°C for 24 hours.

DISCLAIMER

The statements and technical information contained herein are believed to be accurate as of the date hereof. Ossfila Technology Limited reserves the rights to modify, edit, or remove any of said content at our discretion, without notice, and for any reason. Since the conditions and methods of use of the product and of the information referred to herein are beyond our control, Ossfila Technology Limited disclaims any and all liability arising out of the application or use of any product.

SUPPORT

For the latest specifications, additional product information and sales and distribution locations, please contact us at phone: (852) 3188 9234 or email: info@ossfila.com